

DC-7 Diagnostic Ultrasound System

www.mindray.com

MINDRAY



Performance and Value Working for You



MINDRAY™

Mindray Building, Keji 12th Road South, High-tech Industrial Park,
Nanshan, Shenzhen 518057, P. R. China

Tel: +86 755 26582492, 26582888 Fax: +86 755 26582680
E-mail: intl-market@mindray.com Website: www.mindray.com



DISTRIBUTOR:

S S MEDICORP

124, Rajdanga (Gr. Flr), Nabapalli, Kolkata – 107
Phone : +91 33 2441 4244 / 2441 4245 / 4003 1332
Email : infosales@ssmedicorp.com / infoservice@ssmedicorp.com

DC-7 
Diagnostic Ultrasound System



MINDRAY™ is a trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Microstream is a registered trademark or a trademark of ORIDION. MasimoSET is a registered trademark or a trademark of Masimo corporation. Oximax is a registered trademark or trademark of NELLCOR. Specification subject to be changed without notice. ©2009 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.
P/N:ENG-DC-7-420285x4-20090721

DC-7

Diagnostic Ultrasound System

What matters most:

excellent image quality in a value added ultrasound

Based on the latest imaging technologies, the DC-7 is engineered to deliver value added performance. Experience the excellent image quality along with a host of proven capabilities, only seen in premium systems.

- 4D imaging: real time, more details
- Contrast Imaging: enhance pathology structures
- iClear™: speckle noise reduction
- iCape™ view: wide field of view imaging
- iBase inversion tissue harmonic imaging
- iWide bandwidth and high density transducer technology for increased detail and contrast resolution

What is the real ergonomic design: user centric design for optimal operations

multi-function trackball and navigational
primary control are surrounded by 80% of the
systems core operations, minimizing your
fatigue after a long day of scanning.



High resolution LCD monitor

- Fully articulated monitor arm for maximum viewing capability
- iZoom: full screen image display to offer more visibility

Touch screen

- User definable selections and layout
- Multi-language input to meet your specific needs

How to improve your efficiency: streamlined Patient Workflow Circle

Experience the fast response functionality within complex
environment of patient care. The DC-7 is your value and
performance gateway to managing patient healthcare.

Easy to scan

- iTouch™: automatically optimize the image
- iStation™: powerful information management
- Auto calculation: LV trace and IMT
- User-defined packages and exam specific reports

