

Glamatic™

ANAESTHESIA VENTILATOR



SPECIFICATION

- ❖ 7" Touch Screen Display
- ❖ **Modes:** VC, PC, SIMV
- ❖ **Tidal Volume:** 50-1500
- ❖ **Frequency:** 6-40 bpm
- ❖ **IE Ratio:** 4:1 - 1:4
- ❖ **PIP:** 10-60 cm H₂O
- ❖ **ALARM:** High Pressure, Low Pressure and Source Failure.
- ❖ **POWER SUPPLY:** AC 220V

