

**High Frequency Electro Surgical Unit with Separate Bipolar Footswitch :**  
 AYU 400W ALLIGATOR DIGITAL ELECTRO SURGICAL UNIT is High Frequency Electro Surgical Unit Microcontroller based innovative design, together with high reliability makes AYU 400W ALLIGATOR digital Electro Surgical Diathermy unit the most favourable choice for doctors. The technical ability of the unit is combined with patient safety are of prime importance both for the patient and the user. The AYU 400W ALLIGATOR Digital Electro Surgical Diathermy is designed with SMPS power supply provides particular safety and reliability. AYU 400W ALLIGATOR is user friendly electro surgical unit.

**Safety in use :**

The technical conceptions of the AYU 400W ALLIGATOR meet the requirements of highly refined electro surgical techniques all national and international safety standards are confirmed. The micro function allows instant variation in the output power.

**OUTPUT CHARACTERISTICS :**

	MODE	RATED LOAD (OHMS)	MAXIMUM POWER (WATTS)
Monopolar Cut	PURE	300	400
	BLEND-I	300	300
	BLEND-II	300	250
	ENDO CUT	300	200
Monopolar Coag.	DESICCATE	300	150
	FULGURATE	300	250
	FORCE	300	200
	SPRAY	300	120
Bipolar	CUT	100	120
	COAG	100	100
	MACRO	100	70
	MICRO	100	50

Unit Size - LxWxH - 355x330x125 mm  
 Box Size - LxWxH - 430x405x280 mm  
 Net Weight - 5Kg. 700Gm. Gross Weight - 11Kg. 200Gm.

