

# Glamatic™

## ELECTRIC SUCTION MACHINE 1/2 HP



This type of electric suction apparatus adopts completely plastic panel design, which makes it more fancy and fashionable. The device of intercepting oil for circumfluence can reduce the oil pollution. Easy to carry, low noise, high negative pressure and large flux. It can be widely applied in the surgical operations which need to absorb phlegm, thick liquid etc.

- Oil-free pump -Convenient - no need to add oil for maintenance
- Clean - no pollution of oil -Safe - no positive pressure during using

#### TECHNICAL DATA:

Power Voltage: AC220V 22V, 50Hz  $\pm$ 1Hz

Operation noise:  $\leq$  65 dB (A)

Pumping frequency:  $\geq$  20L/min

Max negative pressure:  $\geq$  90kPa

Input power: 180VA

Reservoir capacity: 2500mL/pc, 2pieces

